SCANBODIES

Zfx™ Evolution Matchholder

English

Instructions for Use

DESCRIPTION

Zfx™ Evolution Matchholder (laboratory scan body) is a two-piece scan body with a titanium base and plastic (PEEK) top. The Matchholder is designed to be compatible with a specific implant system and specific diameter. The Matchholder engages the abutment connection of the implant analog and is secured with a captured retaining screw. When scanned, the Matchholder transmits the location and orientation of the implant analog into the design software. The Matchholder features a barcode that corresponds to the implant system and diameter.

INDICATIONS

Zfx™ Evolution Matchholder is indicated for use as a positioning component during the scan of the master cast. Matchholders are indicated for multiple useage.

CONTRAINDICATIONS

Zfx™ Evolution Matchholder is not for clinical use. The scan body should be replaced if there is any noticeable wear, damage or micro-movement when engaging the implant analog, since the condition of the connection is essential to achieve an accurate scan. Do not mark or modify the codes on the Matchholders as this can lead to the marker not being recognized by the scanner software.

LABORATORY PROCEDURE

- Orient the Zfx™ Evolution Matchholder so that the narrow side is to the facial aspect of the master cast.
- Place the Matchholder into the implant analog in the master cast. Rotating the Matchholder slightly will help engage the connection and ensure proper seating of the component.
- If necessary, the soft tissue mask should be removed from the master cast to ensure that the Matchholder is completely seated on the implant analog.
- Once the Matchholder is in place, insert and hand tighten the retaining screw using a Zfx Torque Wrench (Art-No.: ZFX02002067). The Torque Wrench has a clutch that guarantees an equal tightening of each Matchholder. Proceed with the standard scan procedure of the master cast.

GUARANTEE

The mechanical stability of the titanium base and retaining screw is guaranteed as described by law. This guarantee is provided only when all indications and described techniques are strictly followed. These components are subject to continuous development. Therefore, we reserve the right to carry out product modifications without prior notice.

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Distributed by: Zfx GmbH