



A COMPANY OF  
**ZIMMER BIOMET**

# Zfx™ Converter

Get the power of decoding in your own lab

## ZFX™ CONVERTER

The Zfx Converter is a software-module intended for the decoding of 3D information Encode Healing Abutments. The interface is user-friendly and intuitive. The decoded data can be exported as open STL files and afterwards imported into your CAD software. Restorative libraries are included in Zfx CAD software and are also available for Exocad, 3Shape and Dental Wings. Librarians, software-module and the decoding is for free.\*

## ENCODE® IMPRESSION SYSTEM

What is the Encode Impression System:  
A unique 3-in-1 healing abutment, replacing the need for separate impression copings, scan body, and healing abutments – unique codes communicate implant information to facilitate the CAD/CAM design of the final prosthesis – can be scanned with an IOS scanner or traditional impression of the Encode Healing Abutment.

## BENEFITS OF DOWNLOADING

### ZFX CONVERTER

- Decode autonomously
- Immediate results for final prosthesis design
- Increase your business by providing decoding and design services to clinicians who prefer using Encode®

**FREE**  
Software module  
& decoding!\*



**exocad**

**3shape**

**dental wings**



#### FREE DOWNLOAD

Download the Zfx™ Restorative Libraries for free at [www.zfx-dental.com](http://www.zfx-dental.com)


\* Lab can download Zfx Converter for free and free decoding when a GenTek Restoration is used.



# Join the next level of digital impression taking

## INSTALLATION GUIDE

Install Zfx Converter in just a few steps and start decoding Zfx Encode Healing Abutments immediately.



1

### REGISTER & LOGIN

Create a new account or login at  
[my.login.zfx-dental.com/Login](https://my.login.zfx-dental.com/Login)



2

### DOWNLOAD & INSTALL

Download the Zfx Converter  
from the Download Center and  
install the software:  
[my.login.zfx-dental.com/  
Download/Converter](https://my.login.zfx-dental.com/Download/Converter)



3

### START DECODING

Detailed instructions can be found at  
[my.login.zfx-dental.com/  
Help/Converter](https://my.login.zfx-dental.com/Help/Converter)

#### Biomet 3i (UK & Ireland) Ltd.

1 Bell Street, Maidenhead • Berkshire • SL6 1BU  
Tel. +44(0)800 652 1233 • Fax +44(0)800 840 6814  
[www.zimmerbiometdental.com](http://www.zimmerbiometdental.com)

#### Zfx GmbH

Kopernikusstraße 15 • 85221 Dachau, Germany  
Tel. +49 81 31 33 244-0 • Fax +49 81 31 33 244-10  
[www.zfx-dental.com](http://www.zfx-dental.com)