


Zfx Application Manager

Instruction for Use

Version 1.1.12
2014-02-13
Zfx Innovation GmbH



CONTENT

| | |
|-----------------------------|---|
| 1. Login | 2 |
| 2. Home screen | 4 |
| 3. Application | 4 |
| 4. Subscription | 5 |
| 5. Feature Activation | 6 |


Zfx Application Manager

The Zfx Application Manager is a service provided by Zfx Innovation, which is available via PC. This service enables the users to download and administrate software and licences for the Zfx System.

Some applications are free of charge, for others you have to pay for. They can be directly downloaded and installed on the PC.


1. Login

Enter your user name and password for system registration!



In case you're not registered, please click on the registration button.
For direct access, please save the login data.

- 1 The wheel above enables:
- language settings
 - log out (for registration with another account)
 - exit




2. Home screen

The home screen consists of 3 sections, giving a quick information on your state.

3. Application

The section „applications“ informs on:

- 1 the applications, which have already been installed and activated
- 2 the availability of the dongle, connected with the software
- 3 the availability of upgrades, which can be installed.
The synchronization symbol shows the availability of an upgrade, which can be carried-out. Click on the arrow for starting the download.




- 4 Below the application name, the number of the applied version is visible

4. Subscription

The section „subscriptions“ informs on the applications characteristics of the modules and libraries you've acquired.

This section informs on the amount of support hours, which are still available and also on the hours you've already used.



- 1 In this case, the installed version is 5.0.1
- 2 Validity of licence
- 3 Activated CAD Support displaying the remainable support hours
- 4 By clicking on the arrow aside, you get information about the installed application

5. Feature Activation

The section „Feature Activation“ is carried-out by using the Zfx Product Key Card. On this card, you can find the product name including serial number and the according activation code, which needs to be entered into the intended field. The activation code allows a single activation per customer only.

Please enter the activation code of your Zfx product in the textbox below.

2 ——— C7KB8-FW0JI-CRCGY-3S6JQ ———> Activate ✓ ——— 3

You will find the activation code on your Zfx Product Key Card:

1 ———

zimmer dental | Zfx

Please open the Zfx Application Manager and enter the activation code from your PRODUCT KEYCARD (see below) in the "Feature Activation" section. The activation code expires after 1 year. The General Terms and Conditions of Zfx GmbH apply.

Zfx Product Name
ZFX99999999

C7KB8-FW0JI-CRCGY-3S6JQ

Zfx GmbH, Kupernkunstraße 27, 85221 Dachau, Germany
Phone +49 81 31 33 244 - 0, office@zfx-dental.com, www.zfx-dental.com Zahn success formula

- 1 Activation code
- 2 Enter the activation code
- 3 Click on „Activate“